

Phonak Audéo L RIC Order Form

For mild to severe hearing losses

Phonak Audéo L is the comprehensive rechargeable RIC portfolio which includes three designs, four performance levels and five external receivers, including the ActiveVent™ receiver.

Phonak Slim is a design form factor in the Lumity RIC portfolio for mild to severe hearing losses. It is compatible with four external receivers including the ActiveVent™ Receiver.



Step 1: Order details

Account Number :	
Hearing Clinic:	
Delivery Address:	
Contact Name:	Phone Number:
Email Address:	
Date Ordered:	Date Required:

Step 2: Patient information

Name:						
Audiogram: (Required for Intellivent option)						
	kHz	0.25	0.5	1	2	4
	LEFT AC (dB HL)					
	RIGHT AC (dB HL)					

Step 3: Warranty/VIP service

4th year (\$168 per device incl GST) VIP 24 hour service (\$57.75 incl GST) (VIP not available for Titanium products)

Step 4: Choose hearing instrument (see following page for mould order, HSP codes and colour options)

L	R	Model	Part No.	Colour
		L90-R	050-0889	
		L70-R	050-0888	
		L50-R	050-0875	
		L30-R	050-0874	

L	R	Model	Part No.	Colour
		L90-RT	050-1000	
		L70-RT	050-0999	
		L50-RT	050-0998	
		L30-RT	050-1006	

L	R	Model	Part No.	Colour
		L90-RL	050-1002	
		L70-RL	050-1001	
		L50-RL	050-1012	
		L30-RL	050-1011	

	Model	Part No.	Colour
R	Slim L90-R	050-0851	
R	Slim L70-R	050-0849	
L	Slim L90-R	050-0850	
L	Slim L70-R	050-0848	

Step 5(a): Choose receiver option (Sound Delivery System version 4.0 only)

Standard	Medium	Power	UltraPower – cShell only HF3 filter recommended (NOT compatible with Phonak Slim)
S 0-L S 0-R	M 0-L M 0-R	P 0-L P 0-R	UP 0-L UP 0-R
S 1-L S 1-R	M 1-L M 1-R	P 1-L P 1-R	UP 1-L UP 1-R
S 2-L S 2-R	M 2-L M 2-R	P 2-L P 2-R	UP 2-L UP 2-R
S 3-L S 3-R	M 3-L M 3-R	P 3-L P 3-R	UP 3-L UP 3-R

Step 5(b): Choose ActiveVent Receiver option (compatible with Audéo L-R/L-RT/L-RL & Phonak Slim)

Audéo L-RL is compatible with ActiveVent Receiver, however due to the moving parts within ActiveVent and the increased potential for water intrusion into the receiver, it is not recommended for use with Audéo L-RL.

Moderate	(a) Choose earpieces – SlimTip Titanium AV (impression required)	(b) Choose Earpieces – Universal Silicone SlimTip AV (impression not required)
MAV Receiver 0L MAV Receiver 0R	L R	Universal SlimTip AV A-L Universal SlimTip AV C-L
MAV Receiver 1L MAV Receiver 1R		Universal SlimTip AV A-R Universal SlimTip AV C-R
MAV Receiver 2L MAV Receiver 2R	SlimTip Titanium AV	Universal SlimTip AV B-L Universal SlimTip AV D-L
MAV Receiver 3L MAV Receiver 3R		Universal SlimTip AV B-R Universal SlimTip AV D-R

Step 6 (a): Choose dome coupling option (custom ear piece option below)

Dome style	Dome size
Cap Left Right	
Open Left Right	Small Medium Large
Vented Left Right	Small Medium Large
Power Left Right	Small Medium Large

 Cap  Open  Vented  Power

Step 6 (b): Choose custom ear piece option (impression required)

SimTip (removeable earpiece)

L	R
Material	
Hard Acrylic	Hard Acrylic
Soft Silicone	Soft Silicone
Shell colour	
☐ Clear (21)	☐ Clear (21)
Rosa Tint (T)	Rosa Tint (T)
Beige/Pink (26) (hard only)	Beige/Pink (26) (hard only)
Tan (14) (hard only)	Tan (14) (hard only)

cShell Acrylic
(NOT compatible with Phonak Slim)

L	R
Shell colour	
☐ Beige/Pink (26)	☐ Beige/Pink (26)
Clear (21)	Clear (21)
Tan (14)	Tan (14)
Faceplate colour	
☐ Pink (26)	☐ Pink (26)
Cocoa (22)	Cocoa (22)
Brown (28)	Brown (28)
Transparent (13)	Transparent (13)
Tan (14)	Tan (14)
Black Transp. (47)	Black Transp. (47)

cShell Titanium**
(NOT compatible with Phonak Slim)

L	R
Shell colour	
☐ Titanium	☐ Titanium
Faceplate colour	
☐ Black Transp. (47)	☐ Black Transp. (47)

Impressions

New impressions enclosed

Use scanned data
Existing serial number:

SlimTip

cShell Titanium

cShell

Okay to change to cShell if SlimTip
not possible

* Please note cShell Titanium products are not included free of charge with hearing aid order.
S = Standard, will default if no other selection ticked
^ Please note cShell Titanium is only available in Titanium shell colour and Black faceplate colour.

Step 7: Choose venting, wax management and handling options

Venting (mm)*

L	R
AOV Smart vent size calculation. Requires AC thresholds.	
Occluded	
1.0mm	
1.1 - 1.6mm	
1.7 - 2.2mm	
2.3 - 3.0mm	
Open	

*Cerustop default option for SlimTip HA and only option for SlimTip Titanium.

Please note: A coupling code representing actual vent size is printed on custom earpieces. Please enter this code into Phonak Target when programming your client's hearing aids.

Wax management

L	R	SlimTip Hard Acrylic
☐	☐	Cerustop (CS)
		Ext. receiver tube + wax spring (WE)
		None (X)
L	R	SlimTip Soft Silicone
☐	☐	None (X)
L	R	cShell
☐	☐	Cerustop (CS)
		HF3 (recommended for UP Rx) (H3)
		HF4 (H4)
		Ext. receiver tube + wax spring" (WE)

**Not applicable for cShell Titanium.

Handling

L	R
Removal Line (RF)	
Removal Line + 5mm extension (RFE)	
Canal Lock (CL)	
Skeleton Lock (SL)	

Step 8: CROS Fitting?

There are no compatible CROS Transmitters available for the Lumity platform - please consider Paradise platform CROS solutions.

Special requests:

Step 8: Choose accessories

Cleaning kits
Choose one complimentary option

C&C Cleaning Kit (098-0041)
D-Dry+ Electronic Kit (098-0511-00014)

Complimentary charging accessory
(for Audéo L-R/L-RT)

Phonak Charger Ease (075-0080-10)

Complimentary charging accessory
(for Audéo L-RL)

Phonak Charger Case Et Go (075-3017-10)
Phonak Life Charger (075-3023-10)

Complimentary charging accessory
(for Phonak Slim)

Phonak Slim Charger (075-3024-10)

Other accessories

TV Connector (076-3006-06)
Remote Control (076-0065-P5)
PartnerMic (076-4001-P510)
Roger Accessories - please complete the Phonak Roger Order Form

Additional information

HSP codes

Model	HSP codes	Category	Model	HSP codes	Category	Model	HSP codes	Category	Model	HSP codes	Category
L90-R	B766PHO	C2	L90-RT	B770PHO	C2	L90-RL	B762PHO	C2	Phonak Slim L90-R	B783PHO	C2
L70-R	B767PHO	C2	L70-RT	B771PHO	C2	L70-RL	B763PHO	C2	Phonak Slim L70-R	B784PHO	C2
L50-R	B768PHO	C2	L50-RT	B772PHO	C2	L50-RL	B764PHO	C2			
L30-R	B769PHO	C2	L30-RT	B773PHO	C2	L30-RL	B765PHO	C2			

Housing colours

P1 Sand Beige	P3 Sandalwood	P4 Chestnut	P5 Champagne
P6 Silver Gray	P7 Graphite Gray	P8 Velvet Black	H0 Beige

Housing Colours for Phonak Slim only

 2A Silver/black	 3A Graphite/black	 4A Copper/black	 6A Sandalwood/black
---	---	---	---

Audéo L-R/L-RT/L-RL charging accessories

 Phonak Charger Ease 075-0080-10 Compatible with L-R & L-RT	 Phonak Life Charger 075-3023-10 Compatible only with L-RL	 Phonak Charger Case Et Go 075-3017-10 Compatible only with L-RL	 Phonak Slim Charger 075-3024-10 Compatible only with Phonak Slim
---	--	--	---

Models

 Audéo L-R	 Audéo L-RT	 Audéo L-RL	 Slim L-R
---	--	--	--

Other accessories

 TV Connector 076-3006-06	 RemoteControl 076-0065-P5	 PartnerMic™ 076-4001-P510
--	---	--