

# Phonak Serenity Choice Plus Order Form

Phonak Serenity Choice Plus is a range of custom earplugs with different filters for each application.

## Step 1: Order Details

Account Number: \_\_\_\_\_ Contact Name: \_\_\_\_\_ Client Name: \_\_\_\_\_

Hearing Clinic: \_\_\_\_\_ Phone Number: \_\_\_\_\_ Email Address: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

Date Ordered: \_\_\_\_\_ Date Required: \_\_\_\_\_

## Special Instruction

## Step 2: Use case & custom shell preference

Serenity Choice Plus Music					
Material	Style	Shell colour	R	L	Options
Hard Acrylic	Mini				
Soft Silicone	Mini				
Serenity Choice Plus Shooting & Hunting**					
Material	Style	Shell colour	R	L	Options
Hard Acrylic	Mini				
	Full Concha				
Soft Silicone	Mini				
	Full Concha				
Serenity Choice Plus Sleep					
Material	Style	Shell colour	R	L	
Soft Silicone	Mini				
Serenity Choice Plus Work					
Material	Style	Shell colour	R	L	Options
Hard Acrylic	Mini				
	Full Concha				
Soft Silicone	Mini				
	Full Concha				
Serenity Choice Plus Swim					
Material	Style	Shell colour	R	L	Options
Soft Silicone	Full Concha				
Serenity Choice Plus Motorsport					
Material	Style	Shell colour	R	L	Options
Soft Silicone	Mini				
Serenity Choice Plus Fly					
Material	Style	Shell colour	R	L	Options
Hard Acrylic	Mini				
	Full Concha				
Soft Silicone	Mini				
	Full Concha				
Serenity Choice Plus Comfort					
Material	Style	Shell colour	R	L	Options
Soft Silicone	Mini				

## Step 3: Choose filter type


Low		Medium		High	
Filter	Class	Filter	Class	Filter	Class
	2		2		3
	2		2		3
High (impulse)					
Filter					
** Shooting & Hunting is tested to IPIL standards					
No hearing protection					
Filter					
No filter					
Medium		High		Max	
Filter	Class	Filter	Class	Filter	Class
	4		5		5
	4		5		5
	4		4		4
	4		4		4
No hearing protection					
Filter					
		High		Max	
		Filter	Class	Filter	Class
			4		4
Low					
Filter	Class				
	2				
	2				
	2				
	2				
Filter					

## Colours

Soft Silicone



Hard Acrylic



## Style

Mini



Full Concha



## Hearing Protection Classes

Class	L <sub>Aeq8h</sub> (dB)*	Highest Level of noise (dB)*
1	<90	90
2	90 to <95	95
3	90 to <100	100
4	100 to <105	105
5	105 to <110	110

\*From Occupational Noise Exposure - Selection & use of hearing protectors (2003)

## What is too loud?

Excessive exposure to levels above 85 db can cause hearing damage.

140	Gun shot
130	Jackhammer
120	Rock concert
110	Chainsaw
100	Motorsport
90	Hairdryer
80	Plane cabin
70	Medi iPod vol
60	Conversation
50	Air con
40	Rainfall
30	Dishwasher
20	Leaves rustling
10	Whisper
0	Silence