


# Custom Earpiece Order Form

## Receiver-In-Canal Models – Sound Delivery System (SDS) for Lumity, Paradise & Marvel

### Tube Length

- Measure the tube length to select the correct Receiver or SlimTube
- Measure both ears for binaural fitting



### Vent Styles

- We offer two vent solutions for SlimTip (not available for SlimTip Titanium AV):



#### Regular

- Classic vent solution
- Parallel vent (not modifiable in clinic)



#### Cavity


- Optimized for best fit rate
- Easy to enlarge by the fitter

	S Receiver Gain 47(2cc) • MPO 111(2cc)	M Receiver Gain 51(2cc) • MPO 114(2cc)	MAV Receiver Gain 50(2cc) • MPO 119(2cc)	P Receiver Gain 59(2cc) • MPO 122(2cc)	UP Receiver Gain 71(2cc) • MPO 132(2cc)
<b>Custom Earpieces</b>					
cShell Acryl (solid hard)	•	•		•	•
cShell Titanium (solid hard)	•	•		•	•
SlimTip Acryl (hollow hard)	•	•		•	
SlimTip Silicone (solid soft)	•	•		•	
SlimTip Titanium AV (hollow hard)			•		
<b>Instant Fit – Domes</b>					
Cap Dome	•	•		•	
Open Dome	•	•		•	
Vented Dome	•	•		•	
Power Dome	•	•		•	
Universal SlimTip			•		


### Hearing Instrument Compatibility

The SDS earpiece portfolio is compatible with Lumity, Paradise and Marvel hearing aids. It is not backwards compatible to Belong or earlier generation hearing aids. ActiveVent (MAV) is only compatible with Lumity and Paradise rechargeable RICs.


**S Receiver**  
Mild to moderate hearing loss




**M & MAV Receiver**  
Mild to moderately-severe hearing loss



**P Receiver**  
Moderate to severe hearing loss



**UP Receiver**  
Moderately severe to profound hearing loss



cShell Acryl  
(solid hard)



cShell Titanium  
(solid hard)



SlimTip Acryl  
(hollow hard)



SlimTip Silicone  
(solid soft)



SlimTip Titanium AV  
(hollow hard)



ActiveVent  
Receiver (MAV)

# Custom Earpiece Order Form

Date: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
DD MM Year

**Rush Order in 48hr**  
(\$48 charge) N/A for Titanium

PO#: \_\_\_\_\_

Attention to \_\_\_\_\_  
(at Phonak): \_\_\_\_\_

## STEP 1: General Information

Ship to Account Number:

Bill to Account Number:

Phone#:

Company Name:

Address:

Contact Name (Audiologist/Dispenser):

First Name:

Last Name:

Claim Type: AADL BCEHP CNESST Eastern Health

Greenshield ISC (NIHB) Manitoba Health Nisga'a

ODSP REGIE Supp.Health VAC, DND, RCMP

WCB \_\_\_\_\_ province WSIB

Other: \_\_\_\_\_

Claim # (Required): \_\_\_\_\_  
(Please Specify)

L	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
R	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
	250	500	1K	2K	3K	4K

Please include audiogram with all orders. Audiogram required for Acoustically Optimized Venting (AOV).

## STEP 2: Order Information

A) I will be fitting the earpiece(s) to hearing instrument(s):  L Model: \_\_\_\_\_  R Model: \_\_\_\_\_

B) Please send earpiece(s) & hearing instrument(s):  L Model: \_\_\_\_\_ Colour: \_\_\_\_\_  R Model: \_\_\_\_\_ Colour: \_\_\_\_\_

Please also send me the following accessories: \_\_\_\_\_

## STEP 3: Earpiece Order

S = Standard O = Optional

*00 only available in SlimTube	Please indicate Tube length (00*/0/1/2/3)		Receivers					BTE	Wax System					
	L	R	S Receiver	M Receiver	P Receiver	UP Receiver	ActiveVent Receiver <sup>1</sup>	SlimTube	Cerustop	Spring Wax Guard	Receiver Tube (Ext.)	HF3	HF4	No Wax Guard
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Acryl (hollow hard)	L	R	O	S	O			O	S	O	O			O
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Silicone (solid soft)	L	R	O	S	O			O						Required
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Titanium AV (hollow hard)	L	R					Required					Required		
<input type="checkbox"/> L <input type="checkbox"/> R cShell Acryl (solid hard)	L	R	O	S	O	O			S	O	O	O	O	O
<input type="checkbox"/> L <input type="checkbox"/> R cShell Titanium (solid hard)	L	R	O	S	O	O			Required					

All SlimTip receivers, excluding ActiveVent, will be delivered with a Cerustop Wax Guard in the receiver spout. The SlimTip Acryl will be built with Cerustop as an additional wax protection.

ActiveVent will be delivered with an HF3 wax guard in the receiver spout.

<sup>1</sup>ActiveVent Receiver is only available for Lumity & Paradise rechargeable

## STEP 4: Product Options

Faceplate colour options:

- L R Pink
- L R Clear
- L R Cocoa
- L R Brown
- L R Black

Shell colour options:

- L R Pink (Std. on cShell)
- L R Clear (Std. on SlimTip)
- L R Blue (L) & Red (R) (Acryl 4.0 only)
- L R Brown (Silicone only)
- L R Rose Tint (Silicone only)
- L R Black (Sensor cShell only)

cShell Titanium & SlimTip Titanium AV only available in Titanium Gray

Vent options:

- L R Select-a-Vent (SAV)
- L R Open Vent (SlimTip only)
- L R AOV (Agm. required)
- L R Pressure Vent
- L R No Vent\*

\*No Vent is only available for SlimTip Titanium AV

Other options:

- L R Removed Filament
- L R Extended Removal Filament
- L R Soft Coat
- L R Canal Lock\*

\*Same colour as shell

## Additional Information:

---



---



---



---