

Custom Earpiece 4.0 & 6.0 Order Form

Additional Information

STEP 1: General Information

- Please fill in ship to / bill to account information
- Enter patient information
- Select claim type (if applicable) and include claim #

STEP 2: Hearing Aid Information

- Please check box A or B

STEP 3: Earpiece Order

- Select earpiece type
- Indicate tube length (mandatory for cShell)
- Select receiver type
- Select wax system*
- Enter Audiogram (mandatory for AOV)
- Acryl = Acrylic


*All SlimTip receivers, excluding Active Vent, will be delivered with a CeruShield wax guard in the receiver spout. The SlimTip Acryl will be built with Cerustop as an additional wax protection. ActiveVent will be delivered with an HF3 wax guard in the receiver spout

STEP 4: Product Options

- Select colours
- Select vent option
- Select other options


Tube Length

Measure the tube length to select the correct Receiver or SlimTube. Measure both ears for binaural fitting.



Vent Styles

We offer two vent solutions for SlimTip. N/A on SlimTip Titanium AV



Regular:

- Classic vent solution
- Fixed diameter (not enlargable)



Cavity:

- Optimized for best fit rate
- Easy to enlarge by the fitter

Hearing Instrument Compatibility

The SDS 4.0 earpiece portfolio is compatible with Lumity, Paradise and Marvel hearing aids. It is not backwards compatible to Belong or earlier generation hearing aids. Active Vent (MAV) is only compatible with Lumity and Paradise rechargeable RICs. Sensor receiver 6.0 & Sensor cShell are only compatible with Audeo P-R Fit

	Mild to moderate (S)	Mild to moderately-severe (M)	Mild to moderately-severe (Active Vent MAV)	Mild to moderately-severe (M Sensor)	Moderate to severe (P)	Moderate to Severe (P Sensor)	Moderately-severe to profound (UP)
Domes	Cap Dome Open Dome Vented Dome Power Dome	Open Dome Vented Dome Power Dome	Universal SlimTip	Cap Dome Open Dome Clear Vented Dome	Vented Dome Power Dome		
Custom Ear-Pieces	SlimTip 4.0 Acryl (hollow-hard) SlimTip 4.0 Silicone (solid-soft) cShell 4.0 Acryl (solid-hard) cShell 4.0 Titanium (solid-hard)	SlimTip 4.0 Acryl (hollow-hard) SlimTip 4.0 Silicone (solid-soft) cShell 4.0 Acryl (solid-hard) cShell 4.0 Titanium (solid-hard)	SlimTip Titanium AV (hollow, hard)	Sensor cShell 6.0 (solid hard)	SlimTip 4.0 Acryl (hollow-hard) SlimTip 4.0 Silicone (solid-soft) cShell 4.0 Acryl (solid-hard) cShell 4.0 Titanium (solid-hard)	Sensor cShell 6.0 (solid hard)	cShell 4.0 Acryl (solid-hard) cShell 4.0 Titanium (solid-hard)
Fitting Ranges							

