

# Custom Earmold/Swim Plugs Order Form

Additional Information

**PHONAK**  
life is on

## Earmold Pricing (Private Pay)

Earmold ordered with Phonak BTE Hearing Instrument: **\$50**

Earmold only: **\$50** (Plus regular shipping charges)


Swim Plugs: **\$50** (Plus regular shipping charges)

## Earmold Warranty

Phonak earmolds are guaranteed against defects in material or workmanship for one year. Remakes are free within 90 days of purchase. Phonak Swim Plugs have a 90 day warranty. **Earmolds and Swim Plugs cannot be returned for credit. There is no loss and damage warranty on earmolds.**


## Connector QuickSnap

The QuickSnap system is an easy exchangeable earmold tube with connector. Tube exchange was never so easy. No more drilling and glueing, just exchange the tube by a click! Refer to the order form for the tube type and colour options.




## Earmold Areas

Naming of ear areas



## Prelocked Tubing

Standard for all soft earmolds



## Earmold Colours



## Swim Plug Colours



## Default Configuration

- Acrylic Hard Transparent
- 13 Medium Tube
- Medium Canal Length

# Custom Earmold/Swim Plug Order Form

**PHONAK**  
life is on

Date: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
DD MM Year

PO#: \_\_\_\_\_

Attention to  
(at Phonak): \_\_\_\_\_

## Step 1 General Information

Ship to Account Number:

Bill to Account Number:

Phone#:

Company Name:

Address:

Contact Name (Audiologist/Dispenser):

Claim ☐ VAC/RCMP/DND ☐ HCB (NIHB) ☐ REGIE (RAMQ) ☐ CSST ☐ WSIB (ON) ☐ WorkSafe

Type: BC ☐ WCB Alberta ☐ WCB Manitoba ☐ WCB Sask ☐ WCB Yukon ☐ WHSCC

(NL, WSNB, WCB PEI) ☐ Eastern Health ☐ ODDSP ☐ AADL ☐ Other: \_\_\_\_\_

Client Information M ☐ F ☐ ☐ Changes can not be made without calling

User Name: Last \_\_\_\_\_  
(Please print)

First \_\_\_\_\_

☐ Pediatric

Date of Birth: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
DD MM Year

Claim # (Required): \_\_\_\_\_

## Step 2 Order Details

Select whether you wish the BTE(s) to be delivered with the earmold(s) or send earmold(s) only.

☐ Earmold (Private Pay \$50 each) sent with BTE Hearing Instrument

☐ L Model: \_\_\_\_\_  
L Colour Code: \_\_\_\_\_ (eg. P5)

☐ R Model: \_\_\_\_\_  
R Colour Code: \_\_\_\_\_ (eg. P5)

OR

☐ Earmold only (Private Pay \$50 each)  
(Plus regular shipping charges)

Please also send me the following accessories:

OR

☐ Swim Plugs (Private Pay \$50 each)  
(For Swim Plugs please go to Step 5)

## Step 3 Select Earmold Style



Canal Mold  
☐ L ☐ R



Canal Lock  
☐ L ☐ R



Helix Lock  
☐ L ☐ R




Semi Skeleton  
☐ L ☐ R




Skeleton  
☐ L ☐ R



Carved Half Shell  
☐ L ☐ R



Carved Full Shell  
☐ L ☐ R



Full Shell  
☐ L ☐ R

## Step 4 Select Canal Length / Venting

### Canal Length

Left ☐ Short ☒ Med. ☐ Long  
Right ☐ Short ☒ Med. ☐ Long

### Ear Texture

☐ Soft ☐ Avg. ☐ Firm

L						
R						
	250	500	1K	2K	3K	4K

### Speech Scores

☐ dB HL ☐ dB SPL

	Right	Left	Binaural
SRT	_____	_____	_____
MCL	_____	_____	_____
UCL	_____	_____	_____

☐ AOV Acoustically Optimized Vent  
(Audiogram Required)

OR

☐ Vent Size: Customer Specified:  
☐ L \_\_\_\_\_ ☐ R \_\_\_\_\_

## Step 5 Earmold / Swim Plug Materials

### Earmold Materials

	Clear/Transparent (21)	Rose Tint (7)	Brown/Transparent (N)	Black (06)	Blue (07)	Purple (08)	Red (10)	Orange (11)	Green (17)	White (19)	Yellow (20)	Pink (26)	Glossy finish	No Lacquer (clear only)
Acrylic hard (lucite)	<input type="checkbox"/> L <input type="checkbox"/> R	S	O	NA	NA	NA	NA	NA	NA	NA	NA	NA	S	O
Silicone soft 70	<input type="checkbox"/> L <input type="checkbox"/> R	S	O	O	O	O	O	O	O	O	O	O	S	NA
Silicone very soft 40	<input type="checkbox"/> L <input type="checkbox"/> R	S	O	O	O	O	O	O	O	O	O	O	S	NA

### Swim Plug Materials

	Blue (07)	Red (10)	Green (17)	Yellow (20)	Pink (26)	Light Blue (34)
Floatable Silicone	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R

☐ Please indicate if Swim Plug is for a Lyric wearer

## Step 6 Tube / Connector Details

Soft earmolds come standard with prelocked tubing (TRS)

Hard earmolds come standard with glued tubing

Tube/Connectors are only available in clear/transparent



### Tube/Connector details

	00/0/1/2/3	Slim tube length
13 Medium Tube (13M)*	<input type="checkbox"/> L <input type="checkbox"/> R	NA
13 Thick Wall Tube (13T)	<input type="checkbox"/> L <input type="checkbox"/> R	NA
13 Dry Wall Tube (13D)*	<input type="checkbox"/> L <input type="checkbox"/> R	NA
13 Dry Thick Wall Tube (13W)	<input type="checkbox"/> L <input type="checkbox"/> R	NA
SlimTube (STA)	L	R
Power SlimTube II Naïda B & Sky B (SP & UP models only) (PS2)	L	R

\*Not available with Prelocked tubing

## Step 7 Special Instructions