Price & Reference Guide
February 2024



# Table of Content

Phonak Lumity Products	01
Phonak Paradise Products	06
Phonak Marvel Products	. 11
Phonak Belong Products	17
Phonak Vitus Products	21
Accessories	25
Hearing Aid Batteries	29
Online Policy	31

# Lumity Products

# Conversations Shine with Phonak Audéo™ Lumity•



# Phonak Audéo<sup>™</sup> Lumity Feature Summary

# SmartSpeech lechnology

# utoSense OS 5.0

#### 90 Premium

Dynamic Noise Cancellation
Speech Enhancer
SpeechSensor

Ctarra 7- and 2.0

StereoZoom 2.0

**Motion Sensor Hearing** 

UltraZoom

Roger Technology

Comfort in Echo Speech in Car

Speech in Loud Noise

Music

Comfort in Noise

Speech in Noise

Calm Situation

Streaming Speech

Streaming Music

Comfort in Echo (manual)

Speech in Loud Noise

(manual)

Speech in 360° (manual)

**EchoBlock** 

WindBlock

SoundRelax

DuoPhone

Real Ear Sound

SoundRecover2

WhistleBlock

Tinnitus Balance

20 Channels

# Enjoy easier hearing in a variety of situations<sup>1,2,3</sup>with Phonak SmartSpeech Technology in all levels of Lumity

#### 70 Advanced

**Dynamic Noise Cancellation** 

StereoZoom 2.0

**Motion Sensor Hearing** 

UltraZoom

Roger Technology

Comfort in Echo

Music

Comfort in Noise

Speech in Noise

Calm Situation

Streaming Speech

Streaming Music

Speech in Loud Noise (manual)

Speech in 360° (manual)

WindBlock

SoundRelax DuoPhone

Real Ear Sound

SoundRecover2

WhistleBlock

Tinnitus Balance

20 Channels

# 50 Standard

#### **Motion Sensor Hearing**

UltraZoom

Roger Technology

Comfort in Noise

Speech in Noise

Calm Situation

Streaming Speech

Streaming Music

WindBlock

DuoPhone

Real Ear Sound

SoundRecover2

WhistleBlock

Tinnitus Balance

16 Channels

#### 30 Essential

#### UltraZoom

#### Roger Technology

Speech in Noise

Calm Situation

Streaming Speech

Streaming Music

WindBlock

SoundRecover2

WhistleBlock

Tinnitus Balance

12 Channels

#### All levels of Phonak Lumity include:



connectivity





the myPhonak app



Compatibility with Roger and Phonak wireless accessories

Requires ActiveVent Receiver

<sup>\*\*</sup> Compared to a standard acoustic coupling in qui

<sup>\*\*\*</sup> Compared to APD

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG. Inc. and any use of such marks by Sonova At

<sup>1</sup> Appleton, J. (2020) AutoSense OS 4.0 - significantly less listening effort and preferred for speech intelligibility. Phonak Field Study

News, retrieved from www.phonakpro.com/evidence, accessed August 23rd, 2022

Thibodeau L.M. (2020) Benefits in Speech Recognition in Noise with Remote Wireless Microphones in Group Settings. Journal of

<sup>3</sup> Latzel, M. & Hobi S. (2021) Receiver with mechanical vent provides benefit of open and closed acoustics for better speech understanding in noise and naturalness of own voice perception. Phonak Field Study News, retrieved from www.phonakpro.com

# RIC - Phonak Audéo™ Lumity



Rechargeable





Products	Audéo Audéo Audéo L90-R L70-R L50-R			Audéo L30-R	
Receiver	S / M / P / UP				
MRP (₹) kit	7,96,000	4,10,000 2,16,000		1,20,000	
MRP (₹) per unit incl. charger	4,06,500	2,13,500	1,16,500	68,500	









Products	Audéo L90-RL			Audéo L30-RL
Receiver	S / M / P / UP			
MRP (₹) kit	8,32,000	4,44,000	2,51,000	1,54,000
MRP (₹) per unit incl. charger	4,24,500	2,30,500	1,34,000	85,500

# Fitting Range


S-47



#### M-51




#### P-59



Mild to severe hearing loss.






Moderate to profound hearing loss.

# Phonak Audéo<sup>™</sup> L charging options

# **Phonak Charger Ease**

The Phonak Charger Ease is a compact charger with magnets in the charging slots that make it easy to insert and remove Audéo L-R.



MRP (₹) per unit

8,500

#### **Phonak Life Charger**

An inductive charger for Audéo Life devices without a built-in battery.



MRP (₹) per unit

8,500

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Phonak Naida<sup>™</sup> Lumity



# Phonak Naida<sup>™</sup> L-PR


Rechargeable



Products	Naida L90-PR	Naida L70-PR	Naida L50-PR	Naida L30-PR
MRP (₹) kit	7,96,000	4,10,000	2,16,000	1,20,000
MRP (₹) per unit incl. charger	4,04,500	2,11,500	1,14,500	66,500

#### Fitting Range

#### P Receiver-69db



# Phonak Naida $^{TM}$ L-UP



Products	Naida	Naida	Naida	Naida
	L90-UP	L70-UP	L50-UP	L30-UP
MRP (₹) per unit	3,86,000	1,94,000	99,900	52,500

SlimTube 4.0
HE11 680

# Phonak Sky<sup>TM</sup>L




# Phonak Sky<sup>™</sup> L-PR

Rechargeable



Products	Sky L90-PR	Sky L70-PR	Sky L50-PR	Sky L30-PR
MRP (₹) kit	7,96,000	4,10,000	2,16,000	1,20,000
MRP (₹) per unit incl. charger	4,04,500	2,11,500	1,14,500	66,500

#### UP-84db



# Phonak Sky<sup>™</sup> L-UP



Products	Sky	Sky	Sky	Sky	
	L90-UP	L70-UP	L50-UP	L30-UP	
MRP (₹) per unit			99,900	52,500	

# Phonak Charger Combi BTE



- The next generation charger solution for Naída L-PR and Sky L-PR.
- Convenient and easy to use.

6,500

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Paradise Products

# **Phonak Paradise Feature Summary**

## 90 Premium

Comfort in Echo
Speech in Car
Speech in Loud Noise
Music
Comfort in Noise

Speech in Noise

Calm Situation
Streaming Speech

Streaming Music

Comfort in Echo (manual)

Speech in Loud Noise (manual)

Speech in 360° (manual)

UltraZoom Premium

Dynamic Noise Cancellation

Speech Enhancer

EchoBlock

Tap Control\*

StereoZoom

WindBlock

SoundRelax

Motion Sensor Hearing\*

DuoPhone

Real Ear Sound

SoundRecover2

WhistleBlock RogerDirect

Tinnitus Balance

NoiseBlock

20 Channels

70 Advanced

Music
Comfort in Noise

Speech in Noise

**Calm Situation** 

Streaming Speech

Streaming Music

Speech in Loud Noise (manual)

Speech in 360° (manual)

UltraZoom Advanced

EchoBlock
Tap Control\*

StereoZoom

WindBlock

SoundRelax

Motion Sensor Hearing\*

DuoPhone

Real Ear Sound
SoundRecover2

WhistleBlock

RogerDirect

Tinnitus Balance

NoiseBlock

20 Channels

Available in all technology levels: AutoSense OS<sup>TM</sup> 4.0 intelligent automation

# 50 Standard

Comfort in Noise			
Speech in Noise			
Calm Situation			
Streaming Speech			
Streaming Music			
UltraZoom Standard			
Motion Sensor Hearing*			
DuoPhone			
Real Ear Sound			
SoundRecover2			
WhistleBlock			
RogerDirect			
Tinnitus Balance			
NoiseBlock			
16 Channels			

# 30 Essential

Speech i	in Noise
Calm Sit	cuation
Streamin	ng Speech
Streamin	ng Music
UltraZoo	m Essential
SoundRe	cover2
WhistleB	lock
RogerDir	ect
Tinnitus	Balance
NoiseBlo	ock
12 Chann	nels

<sup>\*</sup>Rechargeable models only

# RIC - Phonak Audéo<sup>TM</sup> P



Moderate to profound hearing loss.

Audéo™P-312						Fitting Range
	Products	Audéo P90-312	Audéo P70-312	Audéo P50-312	Audéo P30-312	S Receiver – 47dB
1	Receiver		S / M /	' P / UP		30 40 50 60 770
	MRP (₹) per unit	3,57,000	1,79,000	94,500	50,500	S RECEIVER - 4/dB  als H.  10 10 20 30 40 50 60 70 80 90 100 1125 250 560 1k 2k 4k 8k  Mild to moderately-severe hearing loss.
Audéo™P-R						M receiver –51dB
Rechargeable	Products	Audéo P90-R	Audéo P70-R	Audéo P50-R	Audéo P30-R	all 98.  10 10 10 10 10 10 10 10 10 10 10 11 11
	Receiver		S / M /	' P / UP		90 100
4//	MRP (₹) kit	7,36,000	3,78,000	2,00,000	1,10,000	125 250 500 1k 2k 4k 8k Mild to severe hearing loss.
	MRP (₹) per unit incl. charger	3,74,500	1,95,500	1,06,500	61,500	-
						P receiver - 59dB  ### ### ### ### ### ### ### ### ### #
						UP receiver - 71dB  dB HL  0  10  20  30  40  50  60  70  80  80  80  80  80  80  80  80  8

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Power BTE – Phonak Naida<sup>TM</sup> P




Naida <sup>™</sup> P-PR					
Rechargeable	Products	Naida P90-PR	Naida P70-PR	Naida P50-PR	Naida P30-PR
	MRP (₹) kit	7,36,000	3,78,000	2,00,000	1,10,000
	MRP (₹) per unit incl. charger	3,74,500	1,95,500	1,06,500	61,500

P R	ecei	ver	=	69	dΒ
dB HL					
0					
10					
20					
30					
40					
50					
60					
70					
80			$\neg$		
90				$\overline{}$	
100					4
110					

Fitting Range

Mild	to	moderately-severe hearing	oderately-severe hearing	3
Inss	all	audiometric configurations	udiometric configuration	Ś

Naida <sup>™</sup> P-UP					
	Products	Naida P90-UP	Naida P70-UP	Naida P50-UP	Naida P30-UP
	<b>MRP (₹)</b> per unit	3,57,000	1,79,000	92,500	48,500



<sup>\*</sup>All values are as per 2cc coupler data

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Custom - Phonak Virto<sup>™</sup> P



Virto <sup>™</sup> P-312						
	Products	Virto P90-312	Virto P70-312	Virto P50-312	Virto P30-312	
	Receiver	M / P / SP / UP**				
	MRP (₹) per unit	3,57,000	1,79,000	93,000	49,000	

Virto <sup>™</sup> P-312 NW	0

Products	Virto P90-312 NW 0	Virto P70-312 NW 0	Virto P50-312 NW 0	Virto P30-312 NW 0		
Receiver	M / P / SP					
MRP (₹) per unit	3,57,000	1,79,000	93,000	49,000		

# Virto<sup>™</sup>P-10 NW 0



Products	P90-10 NW 0	P70-10 NW 0	P50-10 NW 0	P30-10 NW 0	
Receiver	M / P / SP				
MRP (₹) per unit	3,57,000	1,79,000	93,000	49,000	


\*\* UP receiver available only in ITC

# Fitting Range




For mild to moderate hearing loss, all audiometric configurations

#### P Receiver = 50dB




For mild to moderatelysevere hearing loss, all audiometric configurations

#### SP Receiver = 60dB



For moderate to severe hearing loss, all audiometric configurations

#### UP Receiver = 70dB



For moderate to profound hearing loss, all audiometric configurations

# Marvel Products

# **Phonak Marvel Feature Summary**

# 90 Premium

30 i icililalii			
Calm Situation			
Speech in Noise			
Comfort in Noise			
Music	70 Advanced		
Speech in Loud Noise	Calm Situation		
Speech in Car	Speech in Noise		
Comfort in Echo	Comfort in Noise		
Streamed speech	Music	50 Standard	
Streamed music	Streamed speech	Calm Situation	
	Streamed music	Speech in Noise	
Streaming programs:		Comfort in Noise	
Roger Direct + mic*	Streaming programs:	Streamed speech	
Partner Mic + mic*	Roger Direct + mic*	Streamed music	30 Essential
Phone cal + mic*	Partner Mic + mic*		Calm Situation
	Phone cal + mic*	Streaming programs:	Speech in Noise
Features:		Roger Direct + mic*	Streamed speech
UltraZoom	Features:	Partner Mic + mic*	Streamed music
SNR-Boost	UltraZoom	Phone cal + mic*	
DuoPhone	SNR-Boost		Streaming programs:
SoundRecover2	DuoPhone	Features:	Roger Direct + mic*
Real Ear Sound	SoundRecover2	UltraZoom	Partner Mic + mic*
WhistleBlock	Real Ear Sound	SNR-Boost	Phone cal + mic*
NoiseBlock	WhistleBlock	DuoPhone	
WindBlock	NoiseBlock	SoundRecover2	Features:
EchoBlock	WindBlock	Real Ear Sound	SoundRecover2
SoundRelax	SoundRelax	WhistleBlock	WhistleBlock
Tinnitus Balance	Tinnitus Balance	Tinnitus Balance	Tinnitus Balance
auto Acclimatization	auto Acclimatization	auto Acclimatization	auto Acclimatization
20 Channels	16 Channels	12 Channels	8 Channels



\*Rechargeable models only

# RIC – Phonak Audéo<sup>TM</sup> M



Audéo <sup>™</sup> M-312				Fitting Range	
	Products	Audéo M50–312	Audéo M30-312	S Receiver = 46dB	
	Receiver	S / M /	P/UP	700 700 700 700 700 700 700 700	
	MRP (₹) per unit	60,000	40,000	For mild to moderately-severe hearing loss	
Audéo™M-R				M Receiver = 50dB    Mathematical Receiver = 50dB   Receiver   10   Receiver	
Rechargeable Rechargeable	Products	Audéo M50-R	Audéo M30-R	For mild to severe hearing loss	
	Receiver	Receiver S / M / P / UP			
	MRP (₹) kit	1,30,000	84,000		
	MRP (₹) per unit incl. charger	71,500	48,500	P Receiver = 58dB	
				For mild to severe hearing loss	
				UP Receiver = 67dB  ### ### ### ### ### ### ### ### ### #	

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Power BTE – Phonak Naida<sup>TM</sup> M



Naida <sup>™</sup> M-SP				Fitting Range
	Products	Naida M50-SP	Naida M30-SP	SP = 81dB
	<b>MRP (₹)</b> per unit	60,000	40,000	30 60 70 80 80 80 80 80 80 80 80 80 80 80 80 80

# BTE - Phonak Bolero<sup>TM</sup> M



Bolero™M				Fitting Range
	Products	Bolero M50-M	Bolero M30-M	M = 63dB
	MRP (₹) per unit	66,500	44,500	20 30 40 50 60 70 80 90 100 110 110 110 110 110 110 110 110
Bolero <sup>™</sup> M-PR				PR = 66dB
Rechargeable	Products	Bolero M50-PR	Bolero M30-PR	10 20 30 40 40 40 40 40 40 40 40 40 40 40 40 40
	MRP (₹) kit	1,30,000	84,000	500 100 110 125 250 500 1k 2k Ak nk Mild to severe hearing loss, all audiometric configurations.
	MRP (₹) per unit incl. charger	71,500	48,500	SlimTube HE HE11 680

<sup>\*1</sup> kit includes 2 hearing aids and a charger

# Pediatric BTE - Phonak Sky<sup>™</sup> M



Sky™M - M*				Fitting Range
	Products	Sky M50-M	Sky M30-M	M = 63dB
	MRP (₹) per unit	66,500	44,500	To T
Sky™M-PR*				
Rechargeable	Products	Sky M50-PR	Sky M30-PR	PR = 66dB
	MRP (₹) kit	1,47,000	1,01,000	10 20 20 20 20 20 20 20 20 20 20 20 20 20
	MRP (₹) per unit incl. charger	80,000	57,000	Mild to severe hearing loss, all audiometric configurations.
CL- TMM CD*				SP = 81dB ಈ ಟ
Sky™M-SP*	Products	Sky M50-SP	Sky M30-SP	69 HL 0 150 250 350 440 550 660 750 600 900 1000
	<b>MRP (₹)</b> per unit	60,000	40,000	Mild to profound hearing loss, all audiometric configurations.
*Available on demand  *1 kit includes 2 hearing aids and	a charger			SlimTube HE HE11 680

# Custom - Phonak Virto<sup>TM</sup> M



Virto <sup>™</sup> M-312					
	Products	Virto M90-312	Virto M70-312	Virto M50-312	Virto M30-312
	Receiver		M / P / S	SP/ UP**	
	MRP (₹) per unit	3,42,000	1,68,000	60,000	40,000

Virto <sup>™</sup> M-312 NW 0	)
-------------------------------	---

Products	Virto M90-312 NW 0	Virto M70-312 NW 0	Virto M50-312 NW 0	Virto M30-312 NW 0
Receiver		M / F	P / SP	
MRP (₹) per unit	3,42,000	1,68,000	60,000	40,000

# Virto<sup>™</sup>M-10 NW 0



Products	Virto M90-10 NW 0	Virto M70-10 NW 0	Virto M50-10 NW 0	Virto M30-10 NW 0
Receiver		M / F	P / SP	
<b>MRP (₹)</b> per unit	3,42,000	1,68,000	60,000	40,000


\*\* UP receiver available only in ITC

# Fitting Range




For mild to moderate hearing loss, all audiometric configurations

#### P Receiver = 50dB




For mild to moderatelysevere hearing loss, all audiometric configurations

#### SP Receiver = 60dB



For moderate to severe hearing loss, all audiometric configurations

#### UP Receiver = 70dB



For moderate to profound hearing loss, all audiometric configurations

# **Belong Products**

# **Phonak Belong Feature Summary**

# 90 Premium

Automatic programs:			
Comfort in Echo			
Speech in Car			
Speech in Loud Noise	70 Advanced		
Music	Automatic programs:		
Comfort in Noise	Music		
Speech in Noise	Comfort in Noise	50 Standard	
Calm Situation	Speech in Noise	Automatic programs:	
	Calm Situation	Comfort in Noise	
Signature features:		Speech in Noise	
UltraZoom - SNR-Boost	Signature features:	Calm Situation	
EchoBlock	UltraZoom - SNR-Boost		
SoundRelax	SoundRelax	Signature features:	
WindBlock	WindBlock	UltraZoom - SNR-Boost	30 Essential
DuoPhone	DuoPhone	DuoPhone	Automatic programs:
FlexControl	FlexControl	FlexControl	Speech in Noise
Flex Volume	Flex Volume	Flex Volume	Calm Situation
Real Ear Sound	Real Ear Sound	Real Ear Sound	
RogerReady	RogerReady	RogerReady	Signature features:
Roger and directional setting	Roger and directional setting	Roger and directional setting	UltraZoom - Essential
SoundRecover2	SoundRecover2	SoundRecover2	SoundRecover2
BroadbandBooster	BroadbandBooster	BroadbandBooster	BroadbandBooster
Tinnitus Balance noise generat	or Tinnitus Balance noise generat	orTinnitus Balance noise generat	orTinnitus Balance noise generator
WhistleBlock	WhistleBlock	WhistleBlock	WhistleBlock
NoiseBlock	NoiseBlock	NoiseBlock	NoiseBlock
QuickSync	QuickSync	QuickSync	QuickSync
User Preference Tuning	User Preference Tuning	User Preference Tuning	User Preference Tuning
auto Acclimatization	auto Acclimatization	auto Acclimatization	auto Acclimatization
20 Channels	16 Channels	12 Channels	8 Channels
Rechargeable models only			

<sup>\*</sup>Rechargeable models only

# Power BTE - Phonak Naida™ B




# Naida<sup>™</sup>B SP / UP



Products	Naida B90	Naida B70	Naida B50	Naida B30
MRP (₹) per unit	2,42,000	1,00,000	44,500	31,000






# Custom - Phonak Virto<sup>™</sup> B



Virto	$^{M}\mathbf{B}$					
A		Products	Virto B90	Virto B70	Virto B50	Virto B30
Virto B-312	Virto B-10	Receiver		M / P / S	SP / UP**	
40		MRP (₹) per unit	2,42,000	1,00,000	44,500	31,000
Virto B-10 O	Virto B-312 NW 0					<u>I</u>





# Pediatric - Phonak Sky<sup>TM</sup> B




# $Sky^{TM}B - M^* / P / SP / UP$

Products	Sky B90	Sky B70	Sky B50	Sky B30
MRP (₹) per unit	2,42,000	1,00,000	44,500	31,000

# M = 68dB P = 71dB SP = 81dB Output Wild to moderately-severe hearing loss, all audiometric configurations Wild to severe hearing loss, all audiometric configurations Wild to severe hearing loss, all audiometric configurations Wild to severe hearing loss, all audiometric configurations

# Sky™ B - RIC

	Products	Sky B90	Sky B70	Sky B50	Sky B30
41)	Receiver	S / P / SP			
	MRP (₹) per unit	2,42,000	1,00,000	45,500	31,500



<sup>\*</sup>For every unilateral fitting, the charges will be charged separately.

# Vitus Products

# **Phonak Vitus Feature Summary**


# Vitus+

Everyday Automatic:	
Calm Situation	
Speech in Noise	
Additional programs:	
Max. additional programs 4	
Speech in Noise	
Calm Situation	Vitus
Comfort in Noise	
Music	Additional programs:
Acoustic phone	Max. additional programs 4
Custom program	Speech in Noise
Wireless Fitting	Calm Situation
	Comfort in Noise
Features:	Acoustic phone
UltraZoom	Custom program
SoundRecover	
User Preference Tuning	Features:
WhistleBlock	UltraZoom
NoiseBlock	User Preference Tuning
QuickSync	WhistleBlock
IP Rating 68	NoiseBlock
Auto Acclimatization	IP Rating 68
6 Channels	4 Channels

# Phonak Vitus+TM



Vitus+ RIC			Fitting Range
- M	Products	Vitus + 312	S Receiver – 56 dB  P Receiver – 65 dB
1	Receiver	S / P / UP	70 70 70 70 70 70 70 70 70 70 70 70 70 7
MRP (₹) per unit	25,500	70 80 70 80 90 90 100 100 100 100 100 100 100 100	
			UP Receiver – 72 dB   dB HL  0  10  20  30  40  50  60  70  80  90  110
			125 280 500 1k 2k 4k 8k  Moderate to severe hearing loss.





Vitus+ ITE				Fitting Range			
		Products	Vitus + 10 NW 0	Vitus + 312	M Receiver = 50	P Receiver = 60	
		Receiver	M / P / SP / UP**		70 70 80 90 70 70	90 00 70 70 80	
		MRP (₹) per unit	25,5	500	For mild to moderate hearing loss, all audiometric configurations	For mild to moderately-severe hearing loss, all audiometric configurations	
					SP Receiver = 70  ### 175	UP Receiver = 78  ### Ht.    0	
		** UP receiver availab	ole only in ITC		hearing loss, all audiometric configurations	hearing loss, all audiometric configurations	

# BTE - Phonak Vitus<sup>TM</sup>



Vitus BTE			Fitting Range	
	Products	Vitus	P = 71dB	UP = 81dB
	Receiver	P / UP	20 30 40 50 60 70 80 90 100 1155 250 500 11 28 44 88  Mild to severe hearing loss, all audiometric configurations	20 30 40 50
	MRP (₹) per unit	19,000		Severe to profound hearing loss, all audiometric configurations.
				SlimTube HE / Power SlimTube Regular tube

Vitus Products 24

# **Phonak Chargers**

# **Phonak Charger Ease**

#### **Compatibility**: Audéo<sup>™</sup> L-R



MRP (₹) per unit

8,500

# Phonak Life Charger



**Compatibility : Audéo**™ L-RL / P-RL

MRP (₹) per unit

8,500

# Phonak Charger Case Combi 2

## Compatibility: Audéo™ P-R / M-R Cros P-R



MRP (₹)
per unit

6,500

# **Phonak Charge and Care**



Compatibility : Audéo™ P-R / M-R

Cros P-R

MRP (₹) per unit

8,500

# Charger Combi BTE 2

 $\begin{array}{c} \textbf{Compatibility: Naida}^{\mathsf{TM}} \ \textbf{L-PR} \ / \ \textbf{P-PR} \\ Sky^{\mathsf{TM}} \ \textbf{L-PR} \ / \ \textbf{M-PR} \end{array}$ 





MRP (₹) per unit

6,500

# Phonak CROS-L



CROS L-R		_		Compatibility			
- A	Products	CROS L-R		Audéo	Na	ida	Sky L-PR
	MRP (₹) per unit	48,500		L-R		PR	
			CROS L-R	•		•	•
CROS P-R				Compatibility			
m)	Products	CROS P -R		Audéo™	Audéo™	Audéo™	Audéo™
4	MRP (₹) per unit	44,500		P-R 90	P-R 70	P-R 50	P-R 30
			CROS P-R	•	•	•	•

# Roger

# Roger On<sup>™</sup> iN - V2

#### (with inbuilt receiver)



Versatile microphone dedicated for all conversation where background noise is present

MRP  $(\mathbf{\xi})$  per unit

90,000

# Roger On - V2



## (without inbuilt receiver)

Versatile microphone dedicated for all conversation where background noise is present

Color: Champagne

MRP (₹) per unit

57,500

# Roger Select<sup>™</sup> iN

#### (with inbuilt receiver)

Adaptive microphone suitable for all listening situations



Color: Champagne

MRP (₹)

**75,000** 

# **Roger Select**

#### (without inbuilt receiver)



Adaptive microphone suitable for all listening situations

Color: Champagne

MRP (₹) per unit

44,000

# Roger receiver

## NeckLoop 03

Roger NeckLoop is a Roger receiver that pairs any hearing device with a T-coil to Roger microphones



Color: Champagne

MRP (₹) per unit

32,000

#### X 03 (Single Receiver)



Roger X is a miniature universal Roger receiver that is compatible with hearing instrument and cochlear implant speech processor.

Color: Champagne

MRP (₹) per unit

16,000

# **Phonak Accessories**

# Phonak Partner Mic<sup>™</sup>



#### Lightweight wireless microphone

**Includes :** USB-C charging cable, US power supply, pouch, lanyard

Color: Champagne (P5)

MRP (₹)

19,000

#### **TV Connector**

#### Enhance the TV listening experience



**Includes :** Micro-USB charging cable, US power supply, Optical cable,

3.5mm audio cable

MRP (₹)

10,000

#### **Phonak Remote Control**



#### Discreet remote control

Includes: Cr2032 battery

Color: Champagne (P5)

MRP (₹) per unit

10,500

# **New Phonak Batteries**

# Improved offering of Zinc-Air Batteries




- Featuring an updated package design, our new Phonak battery blister cards carry a modern look and feel.
- Our HI batteries are mercury-free and adhere to our high quality standards.
- Phonak is a well-established and trusted brand in the hearing care industry.
- Zinc-air batteries are available in the following sizes: 10, 13, 312 and 675.
- Our new Phonak zinc-air batteries offer significantly better battery life and improved performance for wireless streaming and hands-free calling (as compared to current Phonak batteries.)

# Phonak Branded Batteries TFT

(Thin Foil Technology)

Item Description	MRP (₹)
Phonak 675 MF (Pack of 6)	240
Phonak 13 MF (Pack of 6)	240
Phonak 312 MF (Pack of 6)	240
Phonak 10 MF (Pack of 6)	240





world's largest and most modern production plant for environment friendly hearing aid batteries.

# PHONAK Authorised supplier of powerone)) Batteries in India.









Item Description	MRP (₹)
POWERONE 675 MF (Pack of 6)	240
POWERONE 13 MF (Pack of 6)	240
POWERONE 312 MF (Pack of 6)	240
POWERONE 10 MF (Pack of 6)	240

Batteries 30

# Guidelines for online sale and distribution of products of Sonova Hearing India Private Limited

Effective: 1st June 2019

Sonova Hearing India Private Limited ("SIPL"), sells medical devices including hearing aids and wireless hearing accessories (the "Products") in India inter alia which bear the trademark and brand of "Phonak" and "Unitron'.

In the Indian market, SHIPL sells its Products through its authorised distributors and dealers, who provide high standards of best practices, including in hearing assessment, fitting, programming and subsequent rehabilitative services to its patients.

Authorised dealers follow a recommended protocol and have knowledge of SHIPL's technology, skills and / or expertise as regards the suitability, programming and fitting of the Products and software, therefore in an individual appointment they are able to provide patients with numerous value added services including but not limited to, (i.) evaluation of the patient's hearing loss through an audiometric test, (i) select the right technology for the patient as per the patient's needs; (ii) fine-tune the hearing aid fitting with available software based on patient-specific behavioral measures and (iii) provide counseling and aural rehabilitative services relative to the patient's specific needs, including care and maintenance of the hearing aid, (iv.) replacements and (v.) instructions on hearing aid features, such as how to utilize it with wireless systems by applying their knowhow of the Products.

Due to such services and support provided to the patients, it is beneficial to the patients that the Products are provided by only our authorised distributors (who sell to our dealers) and our dealers.

In addition, it is SHIPL's right to protect the sanctity of its distribution channels, including both online and offline sales, in an attempt to safeguard the goodwill and reputation of its Products and its brand image.

Every seller who wishes to sell SH1PL's Products online is required to obtain prior written approval from S1-11PL for the sale of its Products. In this manner, the risks associated with online sales of its Products are minimized by SHIPL.

SH1PL will, without assuming any liability, refuse to accept sales of its Products made through online marketplaces by any seller or by e-commerce platform (not being an authorized distributor or dealer of S1-11131\_, who fall under the subsequent paragraph), and no after sales and warranty may be granted by SH1PL (at its discretion) in relation to such unauthorized sales.

The above mentioned guidelines are also applicable to all authorised distributors and dealers of SHIPL, who are hereby restrained from using, marketing, advertising, displaying or supplying for sale, SHIPL's Products on any online platform without the prior written approval of SHIPL. SHIPL will independently determine noncompliance of the above guidelines. In the event an authorized distributor or a dealer does not comply with the above guidelines, SHIPL may take action against such authorized distributor(s) and dealers, including discontinuation of its engagement with them.


As on 1st June 2019, SHIPL has not authorized any person, proprietorship, firm, company or e-commerce marketplace to sell its Products online.

Anyone found listing, selling, marketing, advertising, displaying or supplying for sale, SHIPL's Products on an any online marketplace or e-commerce platform and displaying SHIPL's images and logos in connection with such sale, without obtaining prior written approval from SHIPL, will be considered as infringing SHIPL's intellectual property rights and may be subject to legal proceedings at its sole risk and cost.

<sup>\*</sup> Price list effective from 1st February 2024

This price list supersedes all previous price lists and is subject to changes without any prior notice.

# EAR Anatomy



# life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we thrive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

www.phonak.com

© 2024 Phonak. All rights reserved.



#### Sonova Hearing India Pvt.Ltd

503-506, A & B Wing, Kanakia Zillion, LBS Marg, BKC Annexe, Near Old Kurla Depo, Kurla (W), Mumbai-400070, India. Phone No.: 022 6782 0666 | Toll Free No.: 1800 120 4401

