

# Phonak Bolero™ B



NEW rechargeable  
hearing aids with  
built-in lithium-ion  
battery

When you can enjoy 24 hours\* of hearing with one simple charge, life is on

A Sonova brand



**PHONAK**  
life is on



# Introducing the new Phonak Bolero B with rechargeable technology

From the world's leading provider of life-changing hearing solutions comes the new Bolero B product line with innovative technologies that make hearing effortless.

## Rechargeable model

- 24 hours\* of hearing with one simple charge
- The quickest charging rechargeable hearing aids ever<sup>1</sup>
- Easy handling, no more hassles of disposable batteries
- Groundbreaking built-in rechargeable lithium-ion battery
- Water and dust resistant

## New AutoSense OS™

- Unmatched hearing performance in everyday listening situations<sup>2</sup>
- Simply adapts to every sound environment automatically – no need to manually adjust the hearing aids<sup>3,4</sup>

\* Expected results when fully charged, and up to 80 minutes wireless streaming time.  
Please refer to [www.phonakpro-us.com/evidence](http://www.phonakpro-us.com/evidence)



# When you can enjoy 24 hours\* of hearing with one simple charge, life is on

Phonak Bolero B-PR features a specifically designed **built-in lithium-ion rechargeable battery** that provides 24 hours\* of hearing with one simple charge.

Phonak Bolero B-PR rechargeable hearing aids combine our latest high-performance technology with **groundbreaking battery technology**.

**No more hassles with disposable batteries** that are small and difficult to handle. Now you can enjoy the freedom and confidence of knowing that Bolero B-PR is conveniently charged and ready to accompany you through the day.

\* Expected results when fully charged, and up to 80 minutes wireless streaming time.  
Please refer to [www.phonakpro-us.com/evidence](http://www.phonakpro-us.com/evidence)

# Smart charging options

Bolero B-PR gives you **24 hours\*** of hearing with a simple **3-hour charge**. If you're in a rush, you can rely on a short 30-minute charge to give you up to 6 hours of full performance. Plus Bolero B-PR has **easy-to-use smart charging options**, to charge up wherever you are, without having to worry about running out of power.

\* Expected results when fully charged, and up to 80 minutes wireless streaming time.  
Please refer to [www.phonakpro-us.com/evidence](http://www.phonakpro-us.com/evidence)



### Phonak Charger Case BTE

- A charger, drying kit and protective hard case all in one
- Includes a cleaning tool



### Phonak Power Pack

- Easily attached to the Phonak Charger Case
- Ideal for short trips where no power source is available



### Phonak Mini Charger BTE

- Most compact charging option



# When your hearing aids adapt to every sound automatically, life is on



Sound is all around and changes constantly depending on where you are. This is what makes listening a real challenge and the task hearing aids have to perform very complex.

What makes Bolero B and Bolero B-PR so special is a **unique operating system** called **AutoSense OS**. It **senses your surroundings** and **adjusts your hearing aids every step of the way**. Just turn them on and the rest is **automatic**.

It doesn't matter where you are or what you're doing, the Bolero B hearing aids were created to simplify your life. And now with the **new rechargeable Bolero B-PR**, you can enjoy **24 hours\*** of hearing with **one simple charge**.

\* Expected results when fully charged, and up to 80 minutes wireless streaming time. Please refer to [www.phonakpro-us.com/evidence](http://www.phonakpro-us.com/evidence)





# More hearing performance in situations where you need it the most

## Understand more in noise

60% improvement in speech understanding by zooming in on a single voice in a noisy environment - versus without hearing aids.<sup>5</sup>

## Better speech understanding

20% better speech understanding thanks to the automatic selection of the best settings for everyday listening situation.<sup>3,4</sup>

## Hearing in car

37% reduction in effort when listening to conversations in car - compared to prior technology.<sup>6</sup>

## Hear quiet voices

10% improvement in soft speech understanding with our newest products.<sup>7</sup>

## Natural music experience

Top rated hearing aid for music sound quality.<sup>8</sup>

Bolero B comes  
in 4 models,  
4 performance levels  
and 9 colors to suit  
your individual needs



Sand Beige  
P1



Sandalwood  
P3



Chestnut  
P4



Champagne  
P5



Silver Gray  
P6



Graphite Gray  
P7



Bolero B-M



Bolero B-P



Bolero B-SP

Includes a groundbreaking rechargeable model



Bolero B-PR  
With T-coil and volume control



Velvet Black  
P8



Alpine White  
T7



Beige  
01



Water and dust resistant IP68

# Phonak wireless accessories to boost your hearing aids

Phonak offers a range of wireless accessories to boost the performance of your hearing aids if needed.

Whether you find yourself having difficulty to understand while watching

TV, talking on the phone or listening to music, we developed wireless solutions to suit your every need. They also allow you to fully participate in noisy situations or when there's distance between you and the speaker.





---

## Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

[www.phonak-us.com](http://www.phonak-us.com)