



**Phonak Vitus™**

The key to confident hearing

A Sonova brand

**PHONAK**  
life is on

# The key to confident hearing

## Confidence in technology

As the world's leading provider in the realm of hearing technology, Phonak is dedicated to creating solutions that are not only innovative but also life-changing.

Offering the largest portfolio of hearing solutions available on the market, we cater to all degrees of hearing loss and create products to suit you. With Phonak Vitus and Vitus+, you have the choice of two quality hearing aid families that are based on technology that really works.

We have combined our years of expertise and innovation into a hearing aid that suits your lifestyle and provides you with hearing confidence wherever you go.

## The choice is yours

Our new portfolio delivers solid performance in a multitude of hearing situations. When you choose Vitus or Vitus+, you can rest assured that you are getting a hearing aid of excellent value combined with:

- Proven Phonak technology
- Quality sound and hearing performance
- Durable – water and dust resistant

### Phonak Vitus

Phonak Vitus provides you with three state-of-the-art, compact Behind-The-Ear models that are ideal for mild to profound hearing loss. With Phonak propriety technology, Vitus delivers when you need it most.


### Phonak Vitus+

Phonak Vitus+ offers a discreet Receiver-In-Canal, three In-The-Ear and three Behind-The-Ear hearing aids with an enhanced feature set. Vitus+ is easy to use with intuitive controls and functions automatically in challenging noisy environments.

Packed in an appealing design, there's a Vitus hearing aid to suit your needs.

- Behind-The-Ear models
- Receiver-In-Canal model
- In-The-Ear models

Your hearing loss and cosmetic preference will determine the most suitable Vitus for you.



## Key features

Proven Phonak technology ensures that you get quality hearing performance. Vitus and Vitus+ are available with some key features to give your hearing confidence a real boost.



### Robust and reliable

Created with robust composite material, Vitus hearing aids are reliable for everyday listening situations. Through rain, snow, heat or dust they are built to meet high quality standards.



### Personalized for your needs

Get your hearing care professional to adjust your hearing aid to best suit your hearing needs.

## Feature comparison

Below is a feature comparison between Vitus and Vitus+ so you and your hearing care professional can find the one that is right for you.

| Feature                                   | Benefit   | Vitus | Vitus+ |
|---|---|-------|--------|
| WhistleBlock                              | No annoying whistling   | •     | •      |
| UltraZoom (Directional microphone option) | Speech intelligibility in noisy environments                          | •     | •      |
| Everyday Automatic                        | Automatic activation of UltraZoom in noisy environments               |       | •      |
| SoundRecover                              | Making high-pitched sounds audible again                              |       | •      |
| QuickSync                                 | One press controls both hearing aids                                  |       | •      |
| Auto Acclimatization                      | Hearing aid automatically adapts over time for comfortable adjustment |       | •      |
| Water resistant IP68 <sup>1,2</sup>       | Water and dust resistant  | •     | •      |
| Roger compatible <sup>2</sup>             | Improve speech understanding over distance and in loud noise          | •     | •      |
| Fine-tuning channels                      | Flexible adjustments to your individual hearing                       | 4     | 6      |
| Color options                             | Matching your preferences   | 2     | 5      |

<sup>1</sup> IP68 rating indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

<sup>2</sup> Except Vitus+ In-The-Ear models

# Life is on

At Phonak, we believe that well-hearing equates to well-being and is essential to living life to the fullest. For more than 70 years, we have remained passionate about creating a world where 'life is on' for everyone. Our innovative hearing solutions are designed for people of all ages and all degrees of hearing loss, to connect socially, thrive mentally and emotionally.

[www.phonak.com](http://www.phonak.com)



028-2223-02/V1.00/2020-12/na © 2020 Sonova AG All rights reserved

