


Your hearing. Your solution.



Name _____

Date _____

Your audiogram



Hearing aid styles

ITE
(In-The-Ear)

RIC
(Receiver-In-Canal)

BTE
(Behind-The-Ear)

Available in a host of hair, skin, fashion and fun colors.

Your hearing needs



☐ One to one



☐ Leisure



☐ TV



☐ Telephone



☐ Classroom /
Conference

Your preferences:

Premium	•	•	•	•	•
Advanced	•	•	•	•	•
Standard	•	•	•	•	
Essential	•	•			

Your recommended performance level:

Boost the performance of your hearing aids with Phonak wireless accessories



**Roger Pen /
Roger Receiver**



**Roger Clip-On Mic /
Roger Receiver**



**Phonak DECT
cordless phone**

Ideal for	<ul style="list-style-type: none"> • Family gatherings • One-to-one conversations • Restaurants and bars • Meetings or lectures 	<ul style="list-style-type: none"> • Family gatherings • One-to-one conversations • Restaurants and bars • Meetings or lectures 	<ul style="list-style-type: none"> • Phone calls
Also helpful for	<ul style="list-style-type: none"> • Listening to music • Watching TV • Watching movies on your computer or tablet • Cell phone calls 	<ul style="list-style-type: none"> • Listening to music • Watching TV • Watching movies on your computer or tablet 	

Your personal guide to better hearing

Notes:



☐ Business meeting



☐ Sports



☐ Restaurant



☐ Social event



☐ Music



☐ Driving



☐ Outdoors

•	•	•	•	•	•	•
•	•	•				



**Phonak TVLink S /
Phonak ComPilot**

- Watching TV



**Phonak ComPilot /
Phonak RemoteMic**

- One-to-one conversations
- Meetings or lectures
- Car



Phonak ComPilot

- Listening to music
- Watching TV
- Watching movies on your computer or tablet
- Cell phone calls



Phonak PilotOne

- Remote control

- Listening to music
- Watching TV
- Watching movies on your computer or tablet
- Cell phone calls
- Remote control

- Listening to music
- Watching TV
- Watching movies on your computer or tablet
- Cell phone calls
- Remote control

- Remote control

Phonak feature overview

		Premium	Advanced	Standard	Essential
Understanding everywhere	 Speech in Wind Allows you to enjoy conversations even on a windy day, for example while playing golf.	•			
	 StereoZoom Helps you focus on a single person speaking directly in front of you in noisy crowd, for example while gathering with your family in a restaurant.	•	•		
	 auto StereoZoom Automatically focuses on the single voice you want to hear from a noisy crowd, for example at a lively social event.	•			
	 ZoomControl Focuses on the speaker you want to hear when you cannot face him/ her, for example in a business meeting around a large table.	•	•		
	 auto ZoomControl Tracks and focuses on the speaker, when you cannot face him/ her, for example while in the car with other passengers.	•			
	 DuoPhone Enjoy hearing the phone conversation in both ears, any time, any where.	•	•	•	
Understanding with comfort	 WhistleBlock Eliminates annoying whistling from your hearing aids.	•	•	•	•
	 WindBlock Suppresses wind noise so you can enjoy spending time outdoors, for example when jogging.	•	•		
	 SoundRelax Cushions sudden loud or unpleasant impulse sounds such as clanging dishes or slamming doors.	•	•		
	 NoiseBlock Reduces annoying background noise like engines or fans without reducing the voices you want to hear.	•	•	•	•
	 EchoBlock Helps you to hear more clearly in rooms or areas where there is an echo, for example in a large public hall.	•			
Understanding everything	 SoundRecover Makes high pitched sounds easier to hear in all situations, for example beginnings and endings of words or birds singing.	•	•	•	•
	 UltraZoom Zooms in on voices of people facing you, while noise from the side and back is reduced. For example during a conversation at a table in a small café with two or three friends.	•	•	•	•
	 SoundFlow Automatically adjusts your hearing aid to suit the environment you are in at that time.	•	•	•	
	 Everyday Automatic Automatically adjusts your hearing aid to provide better speech understanding in noise, for example in a conversation with a friend while shopping.				•
	 QuickSync QuickSync allows you to control both hearing instruments with just one touch of the button.	•	•	•	

Phonak has the right hearing aid for you – regardless of your type of hearing loss, your lifestyle, your personal preferences, your age or your budget.

Premium

Developed to satisfy people, who desire and expect nothing but the best. Thanks to a complete set of trailblazing solutions that uniquely and effectively address all hearing challenges. These hearing aids optimally support understanding in even the most difficult listening situations, for example at a busy social event.

Advanced

The perfect balance between performance and price. An impressive package of features that enhances hearing ability in a variety of listening environments, including those with moderate levels of background noise, for example at a restaurant.

Standard

Ideally suited for people who would like value for money while still expecting state-of-the-art technology. A competitive feature set with the capability to support the hearing demands of daily life with minimal to moderate levels of background noise, for example using the telephone in a public place.

Essential

The Essential products are designed for people looking for an affordable yet effective solution to enhance daily communication abilities while mostly in quiet or minimal background noise, for example one to one conversations.

Choose to enjoy better hearing that matches your lifestyle!

Next appointment

Date

Time

This is your personal record of your solution. Please bring it with you to every appointment. Feel free to add notes of your own.

Your hearing care professional: