

SmartLink+ ZoomLink+ EasyLink+

Dynamic FM for teenagers

Today is a beautiful day



PHONAK

life is on




Your hearing aids might do their job perfectly well when it is quiet, but what about when there is more noise? In these situations Phonak's Dynamic FM systems can help.

When can a Dynamic FM transmitter help?

From a distance

If you are situated a few meters away from a friend, effective communication can be tricky. With a Dynamic FM transmitter however, you can hear the speaker's words clearly in your ear.



On the phone

Your hearing difficulties can make using a cellphone tough, prohibiting you from communicating normally with your peers. SmartLink+ niftily overcomes this problem, connecting wirelessly to any Bluetooth cellphone for effective hands-free use.



Using a computer

In today's multimedia age of MP3s and online video, you expect your media to be fully immersive. With a Dynamic FM transmitter, you can plug into any machine with an audio output, whether PC, Mac, TV or Hi-Fi, to enjoy your audio wirelessly.



SmartLink+

Nothing less than the very best

Designed for those who demand the absolute best, this gadget-friendly Dynamic FM transmitter offers exclusive high-end features, including:

- Bluetooth connectivity for liberating hands-free cellphone use
- Advanced Bluetooth support for effective MP3 playback
- A hearing instrument remote control

The SmartLink+ also features: three intelligent microphone beam settings (*shown below*), Phonak's SoftLanding technology (ensuring that you don't hear a distracting bang when the transmitter is placed on a hard surface), a useful audio input that allows you to plug in and enjoy TV, PC and other media devices wirelessly, and a quick-charging battery.

Microphone settings for SmartLink+ and ZoomLink+



SuperZoom

The optimum setting for very noisy situations such as restaurants. The microphone picks up speech from the front and suppresses noise from all other directions.



Zoom

The standard setting for everyday situations such as face-to-face conversations and meetings.

ZoomLink+

Bringing the sound you want directly into your ears

This versatile Dynamic FM wireless microphone offers you better hearing and understanding across a wide range of situations – from family gatherings and parties to restaurants, shopping malls and lecture halls.

This sleek, attractive device includes Dynamic FM technology for listening in very noisy environments, three intelligent microphone beam settings (*shown below*), Phonak's new protective SoftLanding technology, an audio input for effective wireless TV, PC and audio use, and a quick-charging battery.



Omni

Designed for quieter situations, this 360° setting picks up sound from all around the transmitter.

EasyLink+

Making any conversation easier to follow

Ultra-simple to use and therefore suitable for even the youngest users, EasyLink+ is a one-button Dynamic FM microphone.

With just one press for On/Off switching and a fixed microphone beam setting, EasyLink+ offers the crystal-clear performance benefits of Dynamic FM in addition to Phonak's protective SoftLanding technology, the ability to plug into a TV, PC or other audio device for effective wireless listening, and super-quick battery charging.

Microphone settings for EasyLink+



Zoom

EasyLink+'s default setting is perfect for everyday situations such as face-to-face conversations and meetings. This setting can be re-programmed to SuperZoom or Omni by a hearing care professional.

Welcome to the world of Dynamic FM

Thank you for taking the time to learn more about Phonak's new Dynamic FM transmitters – SmartLink+, ZoomLink+ and EasyLink+.

By using these high-performance hearing products with the appropriate Phonak wireless receiver, you or your teenager will be better able to understand speech in a range of situations, ensuring that today truly is a beautiful day.

The 3 new transmitters introduced here bring Phonak's crystal-clear Dynamic FM sound technology to young adults for the first time, working with FM hearing receivers to give those with hearing, listening, attention and learning-related issues a direct link to the speaker, or anything else they want to hear.



SmartLink+

ZoomLink+

EasyLink+

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

**Interact freely. Communicate with confidence.
Live without limit. Life is on.**