

Naída S V

Believe in the power
to hear more



PHONAK

life is on

Treat yourself to power hearing delight

Naída S V is the new power solution.

Naída S V is a high-performing hearing aid, for people with significant hearing loss. To add to your quality of life by choosing audibility, performance, durability and design, read on to find out more.

More hearing

SoundRecover makes high pitched sounds, which are normally difficult to hear, easier to hear.

Annoying whistling is reduced thanks to powerful WhistleBlock technology.



More freedom no matter where or what

Notable dependability

Naída S V combines performance with great looks and water resistance for dependability.

More connectivity

With Phonak wireless accessories, multimedia entertainment, communication and work are much more fun. TV, telephone, MP3 player, GPS and computers can be enjoyed effortlessly.

Pairing the Phonak TVLink with your Naída S V brings TV sound directly to your ears, so you and your family can adjust the volume separately and enjoy watching TV together.

Best performance for an active life

Your benefits	Features
Understanding speech	
Enjoy a fuller spectrum of audible sounds.	SoundRecover
Understanding in background noise	
Automatically zoom on voices in the front and cut down noise from the sides and the back.	UltraZoom Advanced
Hear well when you can't face the speaker.	ZoomControl
Experience comfortable hearing.	NoiseBlock Advanced
Understanding on the phone and other devices	
Connect directly to mobile phones, TVs, MP3 players, etc.	iCom, TVLink, Click'nTalk
Adapting to changing environments	
Enjoy seamless automatic adaptation.	SoundFlow Advanced
Increase hearing comfort outdoors.	WindBlock
Reduce annoying whistling.	WhistleBlock



Improved hearing in noise

Phonak was first to introduce zoom technology to enhance hearing in challenging hearing situations. Naída S V makes zooming even more effective thanks to further improvements in zoom technology.

Naída S V zoom technology:



UltraZoom Advanced

In noisy situations with several people talking around you, UltraZoom automatically zooms to the voices coming from the front. Noise from the side and back is reduced.



ZoomControl

When you cannot easily face the speaker, e.g. in a car, ZoomControl can be selected to zoom to either side and to the back for improved hearing.

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence. Live without limit. Life is on.

www.phonak-us.com

Ask your hearing care professional:



See what others say:
www.1000reasonsforaida.com

