



DynaMic

The inspiro passaround microphone

Phonak's Dynamic FM technology has changed the face of wireless hearing communications in schools, as the inspiro and MultiTalker Network (MTN) have enabled users to better hear, identify and understand the speech of up to 10 teachers at a time.

Now the DynaMic lets students join in the discussion. DynaMic expands the functionality of inspiro and the MultiTalker Network (MTN) by allowing FM users to hear and understand any speaker that picks up the mic. The result? Everyone in the class can be heard, by everyone.

Designed to work exclusively alongside inspiro and operating as another transmitter on the MultiTalker Network, DynaMic integrates inspiro's high performance transmitter technology into a durable handheld microphone.

DynaMic operates as one or more of 10 transmitters on the MTN (the master device still being inspiro), allowing FM-using students to hear their teachers' words via inspiro and fellow students' comments via the microphone.



Standout DynaMic features:

- **inspiro quality:** employs the same acoustics, electronics, software and digital signal processing technology as Phonak's leading transmitter
- **Child-proof:** DynaMic is easy-to-use and super robust
- **Simple, effective design:** features one On/Mute/Off switch and low-battery warning LED with vibra-alarm
- 21 cm long, with a 15 cm antenna
- Weighs just 110 g / 3.8 oz
- Supplied with a sturdy desk stand

The DynaMic passaround microphone places Phonak's industry-leading transmission technology within arm's reach of every student, ensuring no child is left out of the discussion and confirming Phonak's Dynamic FM as the optimal educational hearing system.

Note: DynaMic operates only with inspiro in MTN. It cannot be used as a standalone transmitter.

www.eSchoolDesk.com
www.phonak.com



Dynamis