

The C&C Line care products

The starter kits for the cleaning and care ...

... of your behind-the-ear hearing aids



... of your behind-the-ear devices with an external receiver,
thin sound tube (SlimTube) or your in-the-ear devices



Your hearing-care professional will be happy to give you detailed suggestions on cleaning your hearing systems.

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely.
Communicate with confidence.
Live without limit.
Life is on.

www.phonak.com

Hearing system care and maintenance

A compact handbook

Art. no. 028-1128-02 Care & Cleaning of your hearing system (2011)



Care of hearing systems

The C&C Line refill packs

for cleaning and drying



Cleansing tablets

The combination of cleansing agents with active oxygen for optimum cleaning of your earmolds. Hygienical and effective.

Contents: 20 pieces



Drying capsules

The reliable way to dry your hearing system. The capsules draw moisture from the hearing aids and the earmolds and bond the particles. Thus the life span of your device is increased considerably.

Contents: 6 pieces

The C&C Line refill packs

for disinfection



Cleansing wipes

Easy and uncomplicated cleaning and disinfection for hearing aids and earmolds. The cleansing wipes are soaked in a water-based, fast-acting cleansing liquid that effectively removes deposits and bacteria.

Contents: 25 pieces



Cleansing stick

The care option with up to 5 days antimicrobial long-term effect for your hearing system. Cleansing stick with water-based, disinfecting cleansing liquid based on nanotechnology - for additional care!

Contents: 15 ml